yofill

Solution Introduction

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The screenshot displays the 'Process Map Editor' interface. At the top, the title bar shows '에너지환경 에너지환경3 | XRD를 이용한 Fe-decorated MnSi1.787Al0.01 소재합성용 Fe-decorated MnSi1.787Al0.01 분말 내 이차상분석' and a search bar with '공정도 검색하기' and a user profile 'yofill'. The left sidebar contains navigation options: '게이트웨이', '공정도', '프로토콜', and '장비템플릿'. Below these are search filters for '재료', '공정', '분석목적', '측정목적', and '장비'. A search bar for '재료' is active, showing a list of items including 'Carbon', 'NCM811', 'NMP', 'PVdF', 'LiNixCoyMnzO2...', and several '12345678' entries. The main workspace features a flowchart with nodes: '재료' (purple) -> '공정' (blue) -> '재료' (purple). From the second '재료' node, two paths emerge: one through '분석목적' (blue) and '측정목적' (green) to '장비' (dark blue) and 'Protocol' (white) boxes. A '정렬' (Sort) menu is open, listing options like '왼쪽 맞춤', '가운데 맞춤', '오른쪽 맞춤', '위쪽 맞춤', '중간 맞춤', '아래쪽 맞춤', '가로 간격 동일', and '세로 간격 동일'. A '문서설정' (Document Settings) panel and a 'Map' mini-map are also visible. The right-hand panel shows details for 'NCM811', including '종류: 재료', '이름: NCM811', '형태: 형태없음', and a table of 'Items' with columns 'Category' and 'Items'. A 'Memo' section is also present.

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The screenshot displays a web-based editor for protocol wiki documents. The interface includes a sidebar on the left with navigation icons and a main editing area on the right. The main area is divided into sections for editing different parts of the document, each with a '변경' (Change) and '삭제' (Delete) button.

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XRD를 이용한 Fe-decorated MnSi_{1.787}AlO_{0.01} 소재합성용 Fe-decorated MnSi_{1.787}AlO_{0.01} 분말 내 이차상분석
Secondary phase analysis in Fe-decorated MnSi_{1.787}AlO_{0.01} powder for synthesis of Fe-decorated MnSi_{1.787}AlO_{0.01} material using XRD

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The screenshot shows the software interface for XRD analysis. The main window title is "XRD를 이용한 Fe-decorated MnSi1.787Al0.01 소재합성용 Fe-decorated MnSi1.787Al0.01 분말 내 이차상분석" (Secondary phase analysis in Fe-decorated MnSi1.787Al0.01 powder for synthesis of Fe-decorated MnSi1.787Al0.01 material using XRD). The left sidebar contains navigation options: 게이트웨이 (Gateway), 공정도 (Process Map), 프로토콜 (Protocol), and 장비템플릿 (Equipment Template). The main area displays a "템플릿 불러오기" (Template Import) dialog box with 9 templates:

- 1. 재료 - 공정 - 재료 (Material - Process - Material)
- 2. 재료 - 공정 - 재료 (Material - Process - Material)
- 3. 재료 - 공정 - 재료 (Material - Process - Material)
- 4. 분석목적 - 장비 - Protocol (Analysis Purpose - Equipment - Protocol)
- 5. 분석목적 - 장비 - Protocol (Analysis Purpose - Equipment - Protocol)
- 6. 분석목적 - 장비 - Protocol (Analysis Purpose - Equipment - Protocol)
- 7. 측정목적 - 장비 - Protocol (Measurement Purpose - Equipment - Protocol)
- 8. 측정목적 - 장비 - Protocol (Measurement Purpose - Equipment - Protocol)
- 9. 측정목적 - 장비 - Protocol (Measurement Purpose - Equipment - Protocol)

At the bottom of the interface, there are buttons for "프로토콜 불러오기" (Load Protocol) and "이슈리스트" (Issue List). The bottom right corner shows "템플릿" (Template) and "위원회 승인요청" (Request for Committee Approval).

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The screenshot displays the '장비템플릿 관리' (Equipment Template Management) interface. On the left, there is a sidebar with navigation options: '게이트웨이', '공정도', '프로토콜', and '장비템플릿'. Below this is a list of equipment types, each with a '템플릿 있음' (Template exists) indicator. The main area shows a tree view of a template structure. The tree is organized into sections: 1.작성자, 2.공정도, 3.용어설명, 4.분석목적 (expanded), 5.시료정보, 6.시료전처리, 7.분석장비, and 8.장비점검. The '4.분석목적' section is expanded, showing sub-items like '분석 목적성', '목적', '배경 설명', '시료 형태', '시료 크기', '시료 양', '기타 특징', '전처리 이름', '전처리 장비명', '전처리 과정', '장비명', '장비설명', '점검명', and '설명'. A tooltip points to a '표면의 결합...' node, stating: '이 분석 프로토콜은 EPMA 장비를 사용하여 Ag-decorated (Bi-...'. On the right, there is a '노드 속성값 설정' (Node Attribute Value Setting) panel with a list of UI components: ParentNode, RadioButton, Checkbox, DirectionInput, Mandatory, ShortValue, LongValue, HiddenNode, and ChildNodeToValue. A '템플릿 저장' (Save Template) button is at the bottom of this panel. The top right corner shows a user profile 'yofill' and a 'Map' dropdown.

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